

“Elastomerics, This is the final s t r e t c h...
Let’s use the best.”



Super Ligature Rings
Ligature Ties
Ligature Bulk Cut
Elastic Placer
Separator Modules
Intra Oral Elastics
Extra Oral Elastics
Latex Free Elastics
Elastomeric Chain
Elastomeric Thread
Elastics Tube
Protective Sleeves

Elastomerics



Super Ligature Rings

Features

- Extended elastic memory
- Improved stain resistance
- High force of restoration
- Various colour including grey, yellow, red, clear
- Ideal ring for holding
- Antibacterial effect due to sliver impregnation
- Packing: 1008 EA/ PKG (42x24)

“Improved module design makes ligating easier. Modules can be pulled off of sticks cleanly without burs or weak spots”.

Ligature Ties


Features

- Available in over 50 vibrant, non-fading colors.
- 42 sticks of 24 modules per package.
- Latex free. Safe for latex sensitive patients.
- One size and color per pack.
- Unmatched stretch & rebound.
- Vibrant colors that resist fading.
- High resistance to staining.
- Packing: 1008 EA/ PKG

GLOW IN THE DARK

Colour
Natural
Violet Glow
Light Green
Pink Glow
Bright Orange

1008 PCS PER PACK




Many more colors available.



MULTI COLOURS

Colour
Grey
Clear
Green
Orange
Pink
Bright Red
Bold Blue


1008 PCS PER PACK



Ligature Bulk Cut

Features

- Absorption free elastomeric resin.
- Can be loaded on canes.
- Very economical.
- Available in two colours, clear & grey.
- Packing 1000 EA/PKG



Elastic Placer



Features


- Simplifies placement and removal of intraoral elastics.
- Simply insert elastic placer and stretch to desired position and release.
- Reverse the procedure for elastic removal.
- Mixed assortment of neon pink, neon blue and neon green.

Separator Modules

Radiopaque

Features

- Loose radiopaque separators eliminate any chance of cross contamination. It will show up on x-ray, eliminating lost separators.
- Individually cut modules packaged loose in a bag.
- Available in Regular & Large Sizes.
- Packing: 1000 EA/ PKT




- Regular Size
○ 0.155" OD x 0.042 thickness
- Large Size
○ 0.193" OD x 0.049 thickness

S-type Separators

Features

- These separators, light blue in color, are rounded on the corner to prevent tears.
- Excellent stretch & rebound makes separation easy.
- Tight contacts are no problem as the separator is easy to insert.
- Soft molded modules do not irritate the tissue.
- 70 modules per ring.
- Packing: 5 Rings EA/ PKT



Intra Oral Elastics

Features

- Color coded strong elastic for intra oral use.
- Packing 1000 EA/PKT .
- Packaged in easy to use zip lock bag.




Actual Sizes

- 3 / 16"
- 5 / 16"
- 1 / 4"
- 5 / 16"L
- 3 / 8"

Extraoral Elastics

Amber coloured, strong elastics for use with face mask.

- Available in natural colour only.
- Packing 1000 EA/PKT .



Latex - Free Elastic

Latex Free Elastics are ideal for avoiding patient sensitivity to latex it is made up of a unique polymers which eliminates the protein and other potential allergy reacting compounds associated with latex.


- Available in natural colour only.




Elastomeric Chain

Features

- Unique heat-curing process enables chain to last longer than others.
- Elasticity remains active up to 6 weeks, assuring faster tooth movement.
- Length -15ft per spool.
- Colors available are clear, grey, tooth color, multicolor & glow.





Spool Dispenser



Elastomeric chain is with excellent stretch and rebound, giving superior elongation for easy application. It is available in three sizes. Closed (no space between rings), Short (short space between rings), Long (medium space between rings).

Elastomeric Thread

Solid elastomeric material in precision extruded thread-form. Exert continues long lasting, predictable force. Moisture, odor & stain resistant. Ideal for rotations, retraction & figure of 8's.



Rough surface ensures that knots will not slip

Elastic Tube

Use when spanning an extraction site on by passing area where teeth are not fully exposed. Also products wire against breakage. Soft, supple material slides easily over archwire without snagging.



Protective Sleeve

Features

- Used to slip over archwire to prevent irritation.
- Grey or clear.
- Odor and stain resistant.
- Fits wires up to .045.

