

“Zeta™ Anchor, the final word in orthodontic anchorage”



Orthodontic Anchorage system


Zeta Anchor™ ZA series
Zeta Anchor™ ZB series



Zeta Anchor™

Orthodontic Mini - Screw NEW


“Zeta Anchor™ features a unique design specially designed for orthodontic use. This anchor is very different from other screws that are made just as a add-on product by dental implant companies. Its a dedicated “Ortho-anchorage mini screw”



Zeta Anchor™ Type A (Series ZA)


Part No.	Diameter	Length
ZA 1406	1.4mm	6mm
ZA 1408	1.4mm	8mm
ZA 1410	1.4mm	10mm
ZA 1606	1.6mm	6mm
ZA 1608	1.6mm	8mm
ZA 1610	1.6mm	10mm

UNIQUE 3 stage screw thread: Zeta Anchor™ has a unique feature in it's 3 stage screw thread. It's a way more superior design as compared to uniform screw thread.



Stage 1

- Sharp cutting edge allows very easy self tapping & screws into bone very easily.



Stage 2

- Gives better grip
- Better adaptability
- Rounded edges



Stage 3

- Wedges and stays firmly in the cortex.
- Space between screws is only 0.4mm.
- This avoids bone absorption due to bacterial infections

Zeta Anchor™ Type B (Series ZB)

Part No.	Diameter	Length
ZB 1406	1.4mm	6mm
ZB 1408	1.4mm	8mm
ZB 1410	1.4mm	10mm
ZB 1606	1.6mm	6mm
Zb1608	1.6mm	8mm
ZB 1610	1.6mm	10mm