

Molar Bands Plain

“Anatomically proportioned for each tooth size”

Features

- Molar bands that anatomically conform to the individual tooth.
- Lingual indent helps for a precise fit and reduces occlusal interference.
- Available in proportionally graduated half sizes for each quadrant.
- Quadrant and sizes are permanently laser-etched on to the band.
- Flattened interproximal that make fitting easier and minimum tooth separation.
- Ideally tempered for easy contouring.
- Permanent laser-etched marking makes identification easy.
- A band design that follows gingival contours, rising on the distal and mesial to avoid soft tissue irritation- offering you fewer trial fittings and less trimming.

*Available for both 6th and 7th molars

ANATOMIC MOLAR BANDS (SIZES)

Small - Medium	32	32+	33	33+	34	34+	35
Medium - Large	35	35+	36	36+	37	37+	38
Large - Extra Large	38+	39	39+	40	40+	41	41+

PREFORMED MOLAR BANDS (PLAIN)


Packing	Description
Individual	Any size or quadrant
Intro - Kit (112 pcs)	2 pcs each of any 14 sizes (x 4 quadrants)
Pro - Kit (280 pcs)	5 pcs each of any 14 sizes (x 4 quadrants)
Master-Kit (560 pcs)	10 pcs each of comprehensive range of 14 sizes (x 4 Quadrants)

NEW


“Super snap™ Molar Bands, micro etched for Maximum retention”



Flat Interproximal Walls




7th molar bands available sizes 31 to 38



Laser Marking

<u>16</u> Tooth	UL Quadrant	35 Perimeter	+ Half Size
--------------------	----------------	-----------------	----------------



Micro-Sandblasted for dependable performance



Easy - to - read laser marking



Bands available in convenient organizer box

To order Kits:
Specify Intro- Kit ,
Pro -Kit or Master Kit.

To order Individual bands:
Specify tooth type ,left or
right, size number and quantity.

Band Coding: Legibly marked on the mesial side:
"U" indicates Upper , "L" indicates Lower, "R" indicates Right,
"L" indicates Left.

Our advantage: Custom ordering is possible!!

All bands, whether plain or prewelded can be ordered
as individuals.

BANDS WITH PREWELDED TUBES

Packing	Description
Individual	Any size or quadrant
Intro Kit (112 pcs)	2 pcs each of any 14 sizes (x 4 quadrants)
Pro - Kit (280 pcs)	5 pcs each of any 14 sizes (x 4 quadrants)


ANATOMIC MOLAR BANDS (SIZES)

Small - Medium	32	32+	33	33+	34	37+	35
Medium - Large	35	35+	36	36+	37	34+	38
Large - Extra Large	38+	39	39+	40	40+	41	41+

*Available for both 6th and 7th molars

NEW

Pre - Welded Band with Button & Cleat



Molar Bands

Pre - Welded

Molar Band With Begg's Tubes


- Sizes from 32+ to 41+
- All Begg's tubes with hook.

Molar Band With Roth Tubes

- Option of upper single, Double or Triple.
- Option of lower single or double..

Molar Band With MBT Tubes

- Option of upper single, Double or Triple.
- Option of lower single or double..



Begg's Tube



Single Tube



Double Tube

Optional Headgear / Axillary Slots



Triple Tube