

krystal ice™ Monocrystalline bracket

Roth & MBT Prescription

“krystal-Ice is a true monocrystalline bracket that combines aesthetics & strength to give the “perfect natural look”



The Clear Choice for Superior performance in aesthetics

Krystal - Ice Clear Brackets are an invaluable product for orthodontists who want an aesthetic solution for their practices, one that will deliver superb results for patients. **In line with this requirement we have launched the truly transparent monocrystalline bracket called krystal-Ice. “Krystal - Ice™ is unmatched in its aesthetic & cutting edge technology.**


Features

- Excellent sliding mechanics - Ultra smooth heat polished surface dramatically reduces friction.
- Maximum tie wing strength and advanced heat-treatment technology means Krystal -Ice brackets are twice as resistant to breakage as compared to other ceramic appliances.
- Truly transparent - These brackets are virtually invisible and works with any colour of teeth.
- Smooth and comfortable - these brackets are Contoured for patient comfort.
- Easy to debond - The ball base technology reduces required debonding force.

Bond Strength

Magnified photographs show the unique base. The krystal-ice™ bonding base is sintered with zirconia valls creating millions of undercuts that mechanically lock with the bracket adhesive. This process allows for secure retention during the length of treatment and provides for predictable debonding at the end of treatment.

Sintered Zirconia balls at the base improves bonding strength



Patient Comfort

Each krystal-ice™ bracket has been diamond and heat polished to produce rounded, smooth corners for uncompromised patient comfort. The small bracket size and low profile design also provide a desirable bracket both clinicians and patients will appreciate. The krystal-ice™ brackets are also functional, the tie wings are easy to ligate, exceptionally strong and highly fracture resistant.

Krystal- ice is available in both MBT & ROTH prescriptions !!

Our advantage: We can give you custom kits !!

Full kits (20 pcs)	$\frac{515}{515}$
Extraction kits	$\frac{414}{414}$
Upper/Lower (3-3)	$\frac{313}{313}$
Upper only (5-5)	$\frac{515}{515}$
Upper (3-3)	$\frac{313}{313}$

Easy Bonding/Debonding



Just squeeze the bracket

