

Universal Buccal Tube


Bondable & Weldable

Features

- Same tube can be used as weldable & bondable
- Low profile design
- Funneled wire guide makes wire entry very easy

Maxillary	Torque	Offset
Double	-14	10°

Mandibular	Torque	Offset
Single	-20	0°
Single	-25	0°



Non-Convertible System Information

Mesial/Distal Dimensions: Single, Double and Triple are 3.7mm
Auxillary slot dimensions .018 main slots have an .018x.031 auxillary
 .022 main slots have an .022x.031 auxillary

These are designed with funneled and reduced entries , on both the archwire slot and the auxillary slot ,resulting in easy archwire insertion which can reduce patient's chair time.

All tube have integral mesiogingival hooks

Micro-etched base gives exceptional bonding even with reduced pad one of mini-tubes. Clearly defined indent allows for "accurate placement"



D-Tech Universal Buccal Tube with bondable & weldable base

Universal Mini Buccal Tube

Maxillary	Torque	Offset
Single	-10	0°
Single	0°	0°

Mandibular	Torque	Offset
Single	-12	0°
Single	-25	0°

Mini System Information

Mesial/Distal Dimensions: 2.41mm
Auxillary slot dimensions .018 main slots have an .018x.031 auxillary
 .022 main slots have an .022x.031 auxillary

